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ANALYSIS OF ASBESTOS IN WATER USING FT-IR AND OPTICAL MICROSCOPY

EU DIRECTIVE 1999/77/EC WAS TRANPOSED INTO SPANISH LAW BY A MINISTERIAL ORDER OF DECEMBER 7 2001. THIS LEGISLATION PROHIBITS THE USE AND SALE OF ALL TYPES OF ASBESTOS AND PRODUCTS CONTAINING ASBESTOS. DESPITE THIS, MANY FACILITIES FEATURE ASBESTOS AND IT IS NECESSARY TO ANALYSE THE FIBRES TO ASSESS ENVIRONMENTAL EXPOSURE TO THIS MATERIAL. IPROMA, A CERTIFIED ANALYSIS LABORATORY, HAS CARRIED OUT A STUDY USING FOURIER-TRANSFORM INFRARED SPECTROSCOPY (FTIR) AND PHASE CONTRAST OPTICAL MICROSCOPY.

Introduction

Asbestos is a group of fibrous metamorphic minerals comprising iron, aluminium, magnesium and calcium silicates. Its different forms belong to the serpentine group (white asbestos) or the amphibole group (blue asbestos), which includes amosite (brown asbestos), antophyllite, tremolite and actinolite fibres, amongst others. Asbestos can be found in the air in the form of simple fibres, fibre aggregates or aggregates of fibres and particles. Observed under the microscope, asbestos fibres appear as flexible, curved or straight threads, depending on their length and nature. The length/width ratio of the fibres tends to be very high and long asbestos fibres have been defined by convention as having a ratio of 3:1, with a maximum diameter of 3µm and a minimum length of 5µm.

Given its insulating, mechanical, chemical, heat and flame resistant properties, added to its relatively low cost, it is understandable that asbestos has numerous industrial applications (products or industrial finishes). Almost all forms of asbestos are highly resistant to acids and alkalis and decompose at high temperatures (800-1,000°C).

Asbestos fibres are released into the air and/or water through the erosion of natural deposits or wear of the products containing them. Small-diameter fibres and small particles that contain fibres can be suspended in the air for long periods and can be carried long distances prior to settling and forming deposits.

The EPA has proposed a maximum concentration of 7·10⁶ fibres/L of water for long fibres (f15µm in length). The US Occupational Safety and Health Administration (OSHA) has established limits of 0.1 fibres/cc. of air for fibres with lengths of f15µm in the workplace (40h/week). In Spain, asbestos is on the list of carcinogenic compounds and is classified in the Limits on

Occupational Exposure to Chemical Agents (V.L.A. – Ambient Limit Values).

EPA Method 100.2 (Determination of Asbestos Structures over 10µm in Length in Drinking Water) is recommended for determining the presence and quantification of asbestos structures in potable water. Asbestos structures are identified

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by morphology, selected area electron diffraction and energy dispersive x-ray analysis by means of transmission electron microscopy. A polycarbonate filter or an MCE membrane filter with a pore size of 0.1µm is used to identify and measure all fibres. For structures >10µm, a pore size of 0.1-0.22µm is sufficient.

The membrane filter/phase contrast optical microscopy method (Spanish National Institute for Safety and Health at Work MTA/MA/051-04 method) applicable to any fibre type (multi-fibre method) is used to determine asbestos fibres in air.

The results obtained represent an index of the numerical concentration of fibres rather than an absolute measure. The application interval depends on acceptable density/concentration requisites of fibres in the samples and the sampling volumes.

The sample is taken by drawing a known volume of air through a membrane filter by means of a sampling pump. The filter is rendered transparent with acetone vapour and is prepared with a liquid contrast medium on a microscope slide. The fibres, which are situated in a certain number of determined fields or areas of the filter, randomly chosen, are counted using a phase contrast microscope (≈500x).

When the fibres have been counted, the number of fibres/sample is calculated using the number of observed fields and the active filter area. With the number of fibres/sample and the volume of air collected, the ambient concentration is obtained (fibras/cc of air).

Iproma analysed the presence or absence of asbestos fibres in a water sample by means of phase contrast optical microscopy, and also by undertaking a quantitative test of asbestos using FTIR, thereby enabling comparison of the results obtained in the two analyses and testing of the methods used.

Methodology

It was first necessary to check the turbidity of the water sample in order to know the degree of obscuration of the water and the suspended solids or colloidal material. Subsequently, two different actions were proposed for the purpose of obtaining the maximum possible information in this respect:

- Detection of the absence/presence of fibres using optical microscopy (counting of fibres).
- Test for the detection of asbestos using FTIR (VDI-3866 method).

Counting of fibres with phase contrast optical microscopy

A fibre count was carried out using phase contrast optical microscopy (Spanish National Institute for Safety and Health at Work MTA/MA/051-04 method).

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Phase contrast optical microscopy is based on the delay that occurs in light waves on passing through objects with different refraction indexes, availing of and amplifying such delays. If the refraction index of the material is greater than that of the surrounding medium, a delay, known as a phase change, occurs. This method induces subtle changes in the refraction index of translucent specimens, enabling very fine details to be seen in the structure. The MTA/MA/051-04 method can give rise to false positives of fibrous objects.

Iproma worked with two aliquots of 10mL and 20mL pre-ozonated sample, prepared in accordance with EPA method 100.2 (Determination of Asbestos Structures over 10µm in Length in Drinking Water), using a reticulated filter with a diameter of 25mm, subsequently rendered transparent with acetone vapour.

A Nikon Eclipse 80i microscope (with phase contrast lenses) was used, verified against a Mark II HSE/NPL phase contrast test slide (500x magnification).

A Walton-Beckett gratiue calibrated against a NIST certified micrometer was used to count the fibres.

Quantification by means of Fourier-Transform Infrared Spectroscopy

FTIR is based on the principles of molecular spectroscopy. The basic principle behind molecular spectroscopy is that the molecules absorb energy from light in specific wavelengths, known as their resonance frequencies (vibration).

An infrared spectrometer works with a small sample placed in an infrared cell, where it is subjected to a source of infrared light. The quantity of light absorbed can be calculated by the difference between the intensity of the light before and after passing through the sample cell.

Iproma used an aliquot of 400mL of pre-ozonated sample, filtered through cellulose acetate filters of 47mm with a pore size of 0.2µm.

The filter was then burnt at 500°C in a muffle furnace to remove organic matter and cellulose, so that the ash contained the inorganic compounds.

The ash was mixed with 300mg of KBr and a lozenge of 13mm was made using a hydraulic press at 10Tm prior to scanning with an FTIR spectrometer.

Figura 3. Espectrofotómetro BIO-RAD–Varian Modelo FTS-2.000 Series. | Figure 3. BIO-RAD–Varian Model FTS-2.000 Series spectrometer

Tabla 2. Resultados para el conteo de fibras mediante microscopía óptica.
Table 2. Results of fibre counting with optical microscopy.

Parámetros Parameters	Resultados Results
Turbidez Turbidity	78 U.N.F.
Fibras contadas Fibres counted	2 y 4 fibras 2 y 4 fibres
Campos contados Fields counted	100 campos/filtros 100 fields/filters
Fibras Fibres	254 y 508 fibras 254 y 508 fibres
Filtros Filters	10 y 20 fibras/c.c. 10 and 20 fibres/c.c.
volumen filtrado Volume filtered	25 c.c. 25 c.c.

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
A BIO-RAD – Varian Model FTS-2.000 Series spectrometer was used and a scan was carried out from 4,000 to 3,500 cm⁻¹ (maximum values for absorption of chrysotile, crocidolite and amosite (most commonly used forms of asbestos)).

Results

The results obtained in the study for the two actions –analysis by phase contrast optical microscopy and FTIR- are presented below.

Gráfica 1. Espectros obtenidos (FT-IR). | Graph 1. Spectrums obtained (FT-IR).

Parametric value limits are not set out in our legislation or by the WHO.

The EPA proposes a maximum of 7,000 fibras/cc for drinking water and the sample analysed complies with this limit.

Quantification by FTIR

The FTIR test, with a detection limit of 0.10mg/L, did not

detect asbestos in the sample analysed, as can be seen from the spectrum presented.

Conclusions

Based on the tests carried out on the wastewater sample with the two techniques implemented at the Iproma laboratory, it can be concluded that:

- The presence of asbestos was not detected by means of FTIR spectroscopy, with a detection limit of 0.10 mg/L
- Fibrous structures, the equivalent of 25 fibres/cc, were detected using optical microscopy, well below the maximum of 7,000 fibres/cc proposed by the EPA.

For the analysis of asbestos in water, it can be concluded that both techniques offer satisfactory results in terms of the objective of the analysis.

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